


medartis®

PRECISION IN FIXATION

PRODUCT INFORMATION

Cannulated MTP Reamer System



APTUS®
Foot

Cannulated MTP* Reamer System

Precise reaming – easy handling

Features and clinical benefits


- Five pairs of reamers to fit almost any MTP-1 joint
- Sharp cutting edges for precise bone shape
- Ring on cone reamers ensures even removal of metatarsal osteophytes
- Perfect add-on to the APTUS Foot portfolio



Clinical Example



Reaming of cone reamer



Reaming of cup reamer



Postoperative X-ray

Clinical case published with the kind permission of Paul Rice, Australia

*Metatarsophalangeal joint


Ordering Information

Cannulated MTP Reamers



A-3030.05

1:2



A-3030.06

1:2

Art. No.	Description	Ø	Shaft End	Pieces/Pkg
A-3030.01	Cup	16 mm	AO Quick Coupling	1
A-3030.03	Cup	18 mm	AO Quick Coupling	1
A-3030.05	Cup	20 mm	AO Quick Coupling	1
A-3030.07	Cup	22 mm	AO Quick Coupling	1
A-3030.09	Cup	24 mm	AO Quick Coupling	1

Art. No.	Description	Ø	Shaft End	Pieces/Pkg
A-3030.02	Cone	16 mm	AO Quick Coupling	1
A-3030.04	Cone	18 mm	AO Quick Coupling	1
A-3030.06	Cone	20 mm	AO Quick Coupling	1
A-3030.08	Cone	22 mm	AO Quick Coupling	1
A-3030.10	Cone	24 mm	AO Quick Coupling	1

Cleaning Stylus, Stainless Steel



1:2

Art. No.	System Size	Ø	Length	Pieces/Pkg
A-8000.12	5.0	1.6 mm	222 mm	1

K-Wires, Stainless Steel




A-5040.41


A-5042.41

Art. No.	Ø	Description	Length	Pieces/Pkg
A-5040.41	1.6 mm	trocár	150 mm	10
A-5042.41	1.6 mm	lancet	150 mm	10


Instrument Case and Tray for MTP Reamers



M-6706



A-6622



A-6623

Art. No.	Description	Size	Pieces/Pkg
M-6706	lid for A-6622	120×240 mm	1
A-6622	instrument case	120×240×37 mm	1
A-6623	instrument tray	113×230×31 mm	1

FOOT-03000006_v2 / © 2017-07, Medartis AG, Switzerland. All technical data subject to alteration.

Caution: The devices are restricted to sale only by or on the order of a physician or a hospital.

MANUFACTURER & HEADQUARTERS

Medartis AG | Hochbergerstrasse 60E | 4057 Basel/Switzerland

P +41 61 633 34 34 | F +41 61 633 34 00 | www.medartis.com

USA

Medartis Inc. | 224 Valley Creek Boulevard, Suite 100 | Exton, PA 19341

P +1 610 961 6101 | Toll free 877 406 BONE (2663) | F +1 610 644 2200

SUBSIDIARIES

Australia | Austria | France | Germany | Mexico | New Zealand | Poland | UK | USA

For detailed information regarding our subsidiaries and distributors, please visit www.medartis.com

