





APTUS® Wrist 2.5

medartis®
PRECISION IN FIXATION


ADAPTIVE II Radius Plates, Volar/Drill Guide Blocks




FPL Radius Plates, Volar/Drill Guide Blocks



Radius Small Fragment Plates




1.5 Radius Hook Plates



Radius Correction Plates, Volar


Plates can also be used for the treatment of fractures




APTUS® Wrist 2.5

medartis®
PRECISION IN FIXATION


Radius Fracture Plates, Volar




Radius Fracture Plates, Extra-Articular, Volar




Radius Frame Plates, Volar



Radius Plates XL, Volar




Radius Plates, Dorsal



APTUS® Wrist 2.5

medartis®
PRECISION IN FIXATION

Distal Ulna Plates




Scale 1:1

Order quantity: 1


Plate thickness (mm)

Wrist Fusion Plates, Dorsal




Scale 1:1

Total Wrist Fusion Plate, Dorsal



Scale 1:1

RSL Fusion Plates, Dorsal



Scale 1:1

Screws, K-Wires, Twist Drills and Countersink

Cortical Screws –
Ø 2.5, HexaDrive 7



Scale 1:1

Art. No.	Length	Pieces/Pkg
A-5700.08	8 mm	5
A-5700.10	10 mm	5
A-5700.11/1	11 mm	1
A-5700.12	12 mm	5
A-5700.13/1	13 mm	1
A-5700.14	14 mm	5
A-5700.15/1	15 mm	1
A-5700.16	16 mm	5
A-5700.18	18 mm	5
A-5700.20	20 mm	5
A-5700.22	22 mm	5
A-5700.24	24 mm	5
A-5700.26	26 mm	5
A-5700.28	28 mm	5
A-5700.30	30 mm	5
A-5700.32	32 mm	5
A-5700.34	34 mm	5

For a package of 1, please add /1 at the end of the article number, for example A-5700.14/1

TriLock Locking Screws –
Ø 2.5, HexaDrive 7



Scale 1:1

Art. No.	Length	Pieces/Pkg
A-5750.08	8 mm	5
A-5750.10	10 mm	5
A-5750.12	12 mm	5
A-5750.14	14 mm	5
A-5750.16	16 mm	5
A-5750.18	18 mm	5
A-5750.20	20 mm	5
A-5750.22	22 mm	5
A-5750.24	24 mm	5
A-5750.26	26 mm	5
A-5750.28	28 mm	5
A-5750.30	30 mm	5
A-5750.32	32 mm	5
A-5750.34	34 mm	5

For a package of 1, please add /1 at the end of the article number, for example A-5750.14/1

TriLock Express Screws –
Ø 2.5, HexaDrive 7



Scale 1:1

Art. No.	Length	Pieces/Pkg
A-5755.14	14 mm	5
A-5755.16	16 mm	5
A-5755.18	18 mm	5
A-5755.20	20 mm	5
A-5755.22	22 mm	5
A-5755.24	24 mm	5

For a package of 1, please add /1 at the end of the article number, for example A-5755.14/1

SpeedTip Screws –
Ø 1.5, HexaDrive 4




Scale 1:1

Art. No.	Length	Pieces/Pkg
A-5210.08	8 mm	5
A-5210.10	10 mm	5
A-5210.12	12 mm	5
A-5210.14	14 mm	5

For a package of 1, please add /1 at the end of the article number, for example A-5210.14/1

K-Wires, Stainless Steel



Scale 1:1

Art. No.	Ø	Description	Length	Pieces/Pkg
A-5040.21	1.2 mm	Trocar	150 mm	10
A-5040.41	1.6 mm	Trocar	150 mm	10
A-5042.21	1.2 mm	Lancet	150 mm	10
A-5042.41	1.6 mm	Lancet	150 mm	10
A-5045.41/1	1.6 mm	Olive	60 mm	1
A-5045.41/4	1.6 mm	Olive	60 mm	4

Twist Drills



Art. No.	Ø	Stop	Length	Shaft End	Pieces/Pkg
A-3713	2.0 mm	40 mm	97 mm	Dental	1
A-3723	2.0 mm	40 mm	97 mm	Stryker J-Latch	1
A-3733	2.0 mm	40 mm	91 mm	AO Quick Coupling	1

Scale 1:1

Twist Drills for Gliding Hole



Art. No.	Ø	Stop	Length	Shaft End	Pieces/Pkg
A-3711	2.6 mm	10 mm	67 mm	Dental	1
A-3721	2.6 mm	10 mm	67 mm	Stryker J-Latch	1
A-3731	2.6 mm	10 mm	61 mm	AO Quick Coupling	1

Scale 1:1

Countersink for Cortical Screws



Art. No.	Ø	Length	Shaft End	Pieces/Pkg
A-3830	3.7 mm	45 mm	AO Quick Coupling	1

Scale 1:1

Instrument Case, Trays and Instruments

A-6602.001

Instrument Tray
APTUS Radius,
Lower

Twist Drill
for Gliding
Holes

A-6602.009
Instrument
Case APTUS
Radius

A-2060
Plate and Screw
Holding Forceps, Angled

A-7009
Hook
«Tönnis»

A-7007
Periosteal Elevator

A-6602.002

Instrument Tray APTUS Radius, Upper

A-7013
Wound Retractor
Mini-Langenbeck

Twist
Drills

A-2013
2.5/2.8 Screwdriver
Blade, HD7, AO

M-6707

Lid for Implant and Instrument Case

A-2721
Drill Guide for
Lag Screws

A-3830
Countersink

A-7012
Bone Holding Forceps

A-7006
Bone Elevator
Mini-Hohmann

A-6010.12/A-6010.16
K-Wire Dispenser Ø 1.2
K-Wire Dispenser Ø 1.6

A-2750
Plate Holding
and Positioning
Instrument

A-2073
Handle
with Quick
Connector, AO

A-2726
Drill Sleeve, Self-Holding

A-2730
Depth Gauge

A-2026
2.5/2.8 Drill Guide,
TriLock^{PLUS}

A-2722
Drill Guide, Scaled

A-2047
2.0–2.8 Plate Bending
Pliers, with Pins

A-2310
1.2/1.5 Screwdriver, HD4,
Self-Holding

→ www.medartis.com/products/aptus/wrist

Medartis, APTUS, MODUS, TriLock, HexaDrive and SpeedTip are registered trademarks of Medartis AG/Medartis Holding AG, 4057 Basel/Switzerland
Caution: The devices are restricted to sale only by or on the order of a physician or a hospital.

CE
0197